

Product Datasheet

Human Signal Recognition Particle 54 kDa Protein ELISA Kit (orb1209162)

Description	Human Signal Recognition Particle 54 kDa Protein, SRP54 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	SRP54
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	55,705
Uniprot ID	P61011
Solubility (25°C)	3200ng/L
Sensitivity	7.87ng/L
Entrez	6729
Expiration Date	Please enquire.