

Product Datasheet

Human Platelet-derived Growth Factor Receptor Like Protein ELISA Kit (orb1209164)

Description	Human Platelet-derived Growth Factor Receptor Like Protein, PDGFRL ELISA Kit
Reactivity	Human
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	PDGFRL
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	41,861
Uniprot ID	Q15198
Solubility (25°C)	32ng/ml
Sensitivity	0.028ng/ml
Entrez	5157
Expiration Date	Please enquire.