

Product Datasheet

Human Phenylalanine-4-hydroxylase ELISA Kit (orb1209166)

Description	Human Phenylalanine-4-hydroxylase, PAH ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	PAH
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	51,862
Uniprot ID	P00439
Solubility (25°C)	2400ng/L
Sensitivity	4.86ng/L
Entrez	5053
Expiration Date	Please enquire.