

Product Datasheet

Human Nadh Dehydrogenase ELISA Kit (orb1209194)

Description	Human Nadh Dehydrogenase, Ubiquinone, Iron-sulfur Protein 8, NDUFS8 ELISA Kit
Reactivity	Human
Range	7-1500ng/L
Tested Applications	ELISA
Target	NDUFS8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,705
Uniprot ID	000217
Solubility (25°C)	1600ng/L
Sensitivity	3.39ng/L
Entrez	4728
Expiration Date	Please enquire.