

Product Datasheet

Human L-amino-acidoxidase ELISA Kit (orb1209199)

Description	Human L-amino-acidoxidase, IL4I1 ELISA Kit
Reactivity	Human
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	IL4I1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	62,881
Uniprot ID	Q96RQ9
Solubility (25°C)	16ng/ml
Sensitivity	0.02ng/ml
Entrez	259307
Expiration Date	Please enquire.