

Product Datasheet

Human UI16 Binding Protein 3 ELISA Kit (orb1209213)

Description	Human UI16 Binding Protein 3, ULBP3 ELISA ELISA Kit
Reactivity	Human
Range	0.5-200ng/ml
Tested Applications	ELISA
Target	ULBP3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,949
Uniprot ID	Q9BZM4
Solubility (25°C)	240ng/ml
Sensitivity	0.31ng/ml
Entrez	79465
Expiration Date	Please enquire.