

Product Datasheet

Human Glycosylphosphatidylinositol-anchored High Density Lipoprotein Binding Protein 1 ELISA Kit (orb1209216)

Description	Human Glycosylphosphatidylinositol-anchored High Density Lipoprotein Binding Protein 1, GPIHBP1 ELIS...
Reactivity	Human
Range	0.5-100ng/ml
Tested Applications	ELISA
Target	GPIHBP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	19,806
Uniprot ID	Q8IV16
Solubility (25°C)	128ng/ml
Sensitivity	0.29ng/ml
Entrez	338328
Expiration Date	Please enquire.