

Product Datasheet

Human Dihydrofolate Reductase ELISA Kit (orb1209223)

Description	Human Dihydrofolate Reductase, DHFR ELISA Kit
Reactivity	Human
Range	0.2-70ng/ml
Tested Applications	ELISA
Target	DHFR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,453
Uniprot ID	P00374
Solubility (25°C)	80ng/ml
Sensitivity	0.092ng/ml
Entrez	1719
Expiration Date	Please enquire.