

Product Datasheet

Human Carcinoembryonic Antigen-related Cell Adhesion Molecule 3 ELISA Kit (orb1209244)

Description	Human Carcinoembryonic Antigen-related Cell Adhesion Molecule 3, CEACAM3 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	CEACAM3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,091
Uniprot ID	P40198
Solubility (25°C)	48ng/ml
Sensitivity	0.048ng/ml
Entrez	1084
Expiration Date	Please enquire.