

Product Datasheet

Human Calcium Calmodulin-dependent Protein Kinase Type 1D ELISA Kit (orb1209257)

Description	Human Calcium Calmodulin-dependent Protein Kinase Type 1D, CAMK1D ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	CAMK1D
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	42,914
Uniprot ID	Q8IU85
Solubility (25°C)	4800ng/L
Sensitivity	9.93ng/L
Entrez	57118
Expiration Date	Please enquire.