

Product Datasheet

Human C4B-binding Protein Beta Chain ELISA Kit (orb1209260)

Description	Human C4B-binding Protein Beta Chain, C4BPB ELISA Kit
Reactivity	Human
Range	1-400ng/mL
Tested Applications	ELISA
Target	C4BPB
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	28,357
Uniprot ID	P20851
Solubility (25°C)	480ng/mL
Sensitivity	0.48ng/mL
Entrez	725
Expiration Date	Please enquire.