

Product Datasheet

Human Bromodomain-containing Protein 4 ELISA Kit (orb1209265)

Description	Human Bromodomain-containing Protein 4, BRD4 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	BRD4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	152,219
Uniprot ID	O60885
Solubility (25°C)	2400ng/L
Sensitivity	4.96ng/L
Entrez	23476
Expiration Date	Please enquire.