

Product Datasheet

Human Beta-galactoside Alpha-2 ELISA Kit (orb1209271)

Description	Human Beta-galactoside Alpha-2, ST6GAL1 ELISA Kit
Reactivity	Human
Range	1-400ng/mL
Tested Applications	ELISA
Target	ST6GAL1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46,605
Uniprot ID	P15907
Solubility (25°C)	480ng/mL
Sensitivity	0.48ng/mL
Entrez	6480
Expiration Date	Please enquire.