

Product Datasheet

Human Beta-defensin 104 ELISA Kit (orb1209275)

Description	Human Beta-defensin 104, DEFB104A ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	DEFB104A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	8,526
Uniprot ID	Q8WTQ1
Solubility (25°C)	4800ng/L
Sensitivity	9.65ng/L
Entrez	503618, 140596
Expiration Date	Please enquire.