

Product Datasheet

Human B-cell Differentiation Antigen CD72 ELISA Kit (orb1209280)

Description	Human B-cell Differentiation Antigen CD72, CD72 ELISA Kit
Reactivity	Human
Range	0.5-100ng/ml
Tested Applications	ELISA
Target	CD72
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	40,220
Uniprot ID	P21854
Solubility (25°C)	128ng/ml
Sensitivity	0.25ng/ml
Entrez	971
Expiration Date	Please enquire.