

Product Datasheet

Human B1 Bradykinin Receptor ELISA Kit (orb1209284)

Description	Human B1 Bradykinin Receptor, BDKRB1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	BDKRB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	40,495
Uniprot ID	P46663
Solubility (25°C)	3200ng/L
Sensitivity	7.62ng/L
Entrez	623
Expiration Date	Please enquire.