

Product Datasheet

Human Cholesterol Side-chain Cleavage Enzyme ELISA Kit (orb1209291)

Description	Human Cholesterol Side-chain Cleavage Enzyme, CYP11A1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	CYP11A1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	60,102
Uniprot ID	P05108
Solubility (25°C)	3200ng/L
Sensitivity	8.07ng/L
Entrez	1583
Expiration Date	Please enquire.