

Product Datasheet

Human Wap Kazal Immunoglobulin Kunitz and Ntr Domain-containing Protein 1 ELISA Kit (orb1209320)

Description	Human Wap Kazal Immunoglobulin Kunitz and Ntr Domain-containing Protein 1, WFIKKN1 ELISA Kit
Reactivity	Human
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	WFIKKN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	58,798
Uniprot ID	Q96NZ8
Solubility (25°C)	48ng/mL
Sensitivity	0.057ng/mL
Entrez	117166
Expiration Date	Please enquire.