

Product Datasheet

Human Vesicle-associated Membrane Protein-associated Protein A ELISA Kit (orb1209326)

Description	Human Vesicle-associated Membrane Protein-associated Protein A, VAPA ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	VAPA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,893
Uniprot ID	Q9P0L0
Solubility (25°C)	4800ng/L
Sensitivity	8.91ng/L
Entrez	9218
Expiration Date	Please enquire.