

Product Datasheet

Human Vesicle-associated Membrane Protein 8 ELISA Kit (orb1209327)

Description	Human Vesicle-associated Membrane Protein 8, VAMP8 ELISA Kit
Reactivity	Human
Range	20-7000ng/L
Tested Applications	ELISA
Target	VAMP8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	11,438
Uniprot ID	Q9BV40
Solubility (25°C)	8000ng/L
Sensitivity	8.96ng/L
Entrez	8673
Expiration Date	Please enquire.