

Product Datasheet

Human Ubiquitin-conjugating Enzyme E2 Variant 2 ELISA Kit (orb1209336)

Description	Human Ubiquitin-conjugating Enzyme E2 Variant 2, UBE2V2 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	UBE2V2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	16,363
Uniprot ID	Q15819
Solubility (25°C)	2400ng/L
Sensitivity	5.23ng/L
Entrez	7336
Expiration Date	Please enquire.