

Product Datasheet

Human Tyrosyl-DNA Phosphodiesterase 1 ELISA Kit (orb1209341)

Description	Human Tyrosyl-DNA Phosphodiesterase 1, TDP1 ELISA Kit
Reactivity	Human
Range	20-7000ng/L
Tested Applications	ELISA
Target	TDP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	68,420
Uniprot ID	Q9NUW8
Solubility (25°C)	8000ng/L
Sensitivity	9.85ng/L
Entrez	55775
Expiration Date	Please enquire.