

Product Datasheet

Human Tyrosine-protein Phosphatase Non-receptor Type Substrate 1 ELISA Kit (orb1209342)

Description	Human Tyrosine-protein Phosphatase Non-receptor Type Substrate 1, SIRPA ELISA Kit
Reactivity	Human
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	SIRPA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	54,967
Uniprot ID	P78324
Solubility (25°C)	24ng/mL
Sensitivity	0.031ng/mL
Entrez	140885
Expiration Date	Please enquire.