

Product Datasheet

Human Tyrosine-protein Kinase Receptor Tyro3 ELISA Kit (orb1209346)

Description	Human Tyrosine-protein Kinase Receptor Tyro3, TYRO3 ELISA Kit
Reactivity	Human
Range	0.2-60ng/mL
Tested Applications	ELISA
Target	TYRO3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	96,905
Uniprot ID	Q06418
Solubility (25°C)	64ng/mL
Sensitivity	0.121ng/mL
Entrez	7301
Expiration Date	Please enquire.