

Product Datasheet

Human Tripartite Motif-containing Protein 72 ELISA Kit (orb1209361)

Description	Human Tripartite Motif-containing Protein 72, TRIM72 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	TRIM72
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	52,731
Uniprot ID	Q6ZMU5
Solubility (25°C)	3200ng/L
Sensitivity	7.31ng/L
Entrez	493829
Expiration Date	Please enquire.