

Product Datasheet

Human T-complex Protein 1 Subunit Zeta ELISA Kit (orb1209391)

Description	Human T-complex Protein 1 Subunit Zeta, CCT6A ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	CCT6A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	58,024
Uniprot ID	P40227
Solubility (25°C)	3200ng/L
Sensitivity	7.65ng/L
Entrez	908
Expiration Date	Please enquire.