

Product Datasheet

Human T-cell Surface Protein Tactile ELISA Kit (orb1209394)

Description	Human T-cell Surface Protein Tactile, CD96 ELISA Kit
Reactivity	Human
Range	0.05-18ng/mL
Tested Applications	ELISA
Target	CD96
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	65,634
Uniprot ID	P40200
Solubility (25°C)	19.2ng/mL
Sensitivity	0.027ng/mL
Entrez	10225
Expiration Date	Please enquire.