

Product Datasheet

Human T-cell Immunoglobulin and Mucin Domain-containing Protein 4 ELISA Kit (orb1209396)

Description	Human T-cell Immunoglobulin and Mucin Domain-containing Protein 4, TIMD4 ELISA Kit
Reactivity	Human
Range	0.05-18ng/mL
Tested Applications	ELISA
Target	TIMD4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	41,578
Uniprot ID	Q96H15
Solubility (25°C)	19.2ng/mL
Sensitivity	0.031ng/mL
Entrez	91937
Expiration Date	Please enquire.