

Product Datasheet

Human Succinate-semialdehyde Dehydrogenase ELISA Kit (orb1209409)

Description	Human Succinate-semialdehyde Dehydrogenase, ALDH5A1 ELISA Kit
Reactivity	Human
Range	20-3800ng/L
Tested Applications	ELISA
Target	ALDH5A1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	57,215
Uniprot ID	P51649
Solubility (25°C)	4000ng/L
Sensitivity	9.39ng/L
Entrez	7915
Expiration Date	Please enquire.