

Product Datasheet

Human Steroid Hormone Receptor Err1 ELISA Kit (orb1209412)

Description	Human Steroid Hormone Receptor Err1, ESRRA ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	ESRRA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	45,510
Uniprot ID	P11474
Solubility (25°C)	4800ng/L
Sensitivity	9.84ng/L
Entrez	2101
Expiration Date	Please enquire.