

Product Datasheet

Human Sperm Surface Protein Sp17 ELISA Kit (orb1209420)

Description	Human Sperm Surface Protein Sp17, SPA17 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	SPA17
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	17,406
Uniprot ID	Q15506
Solubility (25°C)	4800ng/L
Sensitivity	10.16ng/L
Entrez	53340
Expiration Date	Please enquire.