

Product Datasheet

Human Sodium-dependent Serotonin Transporter ELISA Kit (orb1209428)

Description	Human Sodium-dependent Serotonin Transporter, SLC6A4 ELISA Kit
Reactivity	Human
Range	50-10000ng/L
Tested Applications	ELISA
Target	SLC6A4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	70,325
Uniprot ID	P31645
Solubility (25°C)	12800ng/L
Sensitivity	30.44ng/L
Entrez	6532
Expiration Date	Please enquire.