

Product Datasheet

Human Smad Nuclear-interacting Protein 1 ELISA Kit (orb1209434)

Description	Human Smad Nuclear-interacting Protein 1, SNIP1 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	SNIP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	45,778
Uniprot ID	Q8TAD8
Solubility (25°C)	4800ng/L
Sensitivity	9.78ng/L
Entrez	79753
Expiration Date	Please enquire.