

Product Datasheet

Human Serine Protease Inhibitor Kazal-type 7 ELISA Kit (orb1209460)

Description	Human Serine Protease Inhibitor Kazal-type 7, SPINK7 ELISA Kit
Reactivity	Human
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	SPINK7
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	9,232
Uniprot ID	P58062
Solubility (25°C)	24ng/mL
Sensitivity	0.029ng/mL
Entrez	84651
Expiration Date	Please enquire.