

Product Datasheet

Human Secreted Ly-6 Upar-related Protein 1 ELISA Kit (orb1209461)

Description	Human Secreted Ly-6 Upar-related Protein 1, SLURP1 ELISA Kit
Reactivity	Human
Range	2-600ng/ml
Tested Applications	ELISA
Target	SLURP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	11,186
Uniprot ID	P55000
Solubility (25°C)	640ng/ml
Sensitivity	1.08ng/ml
Entrez	57152
Expiration Date	Please enquire.