

Product Datasheet

Human Ribonucleoside-diphosphate Reductase Large Subunit ELISA Kit (orb1209471)

Description	Human Ribonucleoside-diphosphate Reductase Large Subunit, RRM1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	RRM1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	90,070
Uniprot ID	P23921
Solubility (25°C)	3200ng/L
Sensitivity	7.58ng/L
Entrez	6240
Expiration Date	Please enquire.