

Product Datasheet

Human Repetin ELISA Kit (orb1209476)

Description	Human Repetin, RPTN ELISA Kit
Reactivity	Human
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	RPTN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	90,731
Uniprot ID	Q6XPR3
Solubility (25°C)	96ng/mL
Sensitivity	0.149ng/mL
Entrez	126638
Expiration Date	Please enquire.