

Product Datasheet

Human Receptor-type Tyrosine-protein Phosphatase Zeta ELISA Kit (orb1209478)

Description	Human Receptor-type Tyrosine-protein Phosphatase Zeta, PTPRZ1 ELISA Kit
Reactivity	Human
Range	5-1000ng/L
Tested Applications	ELISA
Target	PTPRZ1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	254,587
Uniprot ID	P23471
Solubility (25°C)	12800ng/L
Sensitivity	2.48ng/L
Entrez	5803
Expiration Date	Please enquire.