

Product Datasheet

Human Receptor-type Tyrosine-protein Phosphatase F ELISA Kit (orb1209480)

Description	Human Receptor-type Tyrosine-protein Phosphatase F, PTPRF ELISA Kit
Reactivity	Human
Range	2-700ng/L
Tested Applications	ELISA
Target	PTPRF
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	212,879
Uniprot ID	P10586
Solubility (25°C)	800ng/L
Sensitivity	0.83ng/L
Entrez	5792
Expiration Date	Please enquire.