

Product Datasheet

Human Ras-related Protein Rab-27B ELISA Kit (orb1209485)

Description	Human Ras-related Protein Rab-27B, RAB27B ELISA Kit
Reactivity	Human
Range	20-3500ng/L
Tested Applications	ELISA
Target	RAB27B
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	24,608
Uniprot ID	000194
Solubility (25°C)	4000ng/L
Sensitivity	9.04ng/L
Entrez	5874
Expiration Date	Please enquire.