

Product Datasheet

Human Ras-related Protein Rab-14 ELISA Kit (orb1209486)

Description	Human Ras-related Protein Rab-14, RAB14 ELISA Kit
Reactivity	Human
Range	20-3500ng/L
Tested Applications	ELISA
Target	RAB14
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,897
Uniprot ID	P61106
Solubility (25°C)	4000ng/L
Sensitivity	9.81ng/L
Entrez	51552
Expiration Date	Please enquire.