

Product Datasheet

Human Pyrin ELISA Kit (orb1209492)

Description	Human Pyrin, MEFV ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	MEFV
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	86,444
Uniprot ID	O15553
Solubility (25°C)	3200ng/L
Sensitivity	7.96ng/L
Entrez	4210
Expiration Date	Please enquire.