

Product Datasheet

Human Protocadherin-10 ELISA Kit (orb1209505)

Description	Human Protocadherin-10, PCDH10 ELISA Kit
Reactivity	Human
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	PCDH10
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	112,936
Uniprot ID	Q9P2E7
Solubility (25°C)	48ng/mL
Sensitivity	0.056ng/mL
Entrez	57575
Expiration Date	Please enquire.