

Product Datasheet

Human Protein Wnt-7B ELISA Kit (orb1209509)

Description	Human Protein Wnt-7B, WNT7B ELISA Kit
Reactivity	Human
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	WNT7B
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,327
Uniprot ID	P56706
Solubility (25°C)	24ng/mL
Sensitivity	0.022ng/mL
Entrez	7477
Expiration Date	Please enquire.