

Product Datasheet

Human Protein S100-A16 ELISA Kit (orb1209515)

Description	Human Protein S100-A16, S100A16 ELISA Kit
Reactivity	Human
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	S100A16
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	11,801
Uniprot ID	Q96FQ6
Solubility (25°C)	24ng/mL
Sensitivity	0.019ng/mL
Entrez	140576
Expiration Date	Please enquire.