

Product Datasheet

Human Protein Hexim1 ELISA Kit (orb1209521)

Description	Human Protein Hexim1, HEXIM1 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	HEXIM1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	40,623
Uniprot ID	094992
Solubility (25°C)	4800ng/L
Sensitivity	9.61ng/L
Entrez	10614
Expiration Date	Please enquire.