


Product Datasheet

Human Protein Fam13A ELISA Kit (orb1209526)

Description	Human Protein Fam13A, FAM13A ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	FAM13A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	116,932
Uniprot ID	O94988
Solubility (25°C)	3200ng/L
Sensitivity	6.99ng/L
Entrez	10144
Expiration Date	Please enquire.



Standard Curve of Human Protein Fam13A E...