

Product Datasheet

Human Protein Delta Homolog 1 ELISA Kit (orb1209530)

Description	Human Protein Delta Homolog 1, DLK1 ELISA Kit
Reactivity	Human
Range	20-6000ng/L
Tested Applications	ELISA
Target	DLK1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	41,300
Uniprot ID	P80370
Solubility (25°C)	6400ng/L
Sensitivity	9.47ng/L
Entrez	8788
Expiration Date	Please enquire.