

Product Datasheet

Human Protein Cuta ELISA Kit (orb1209532)

Description	Human Protein Cuta, CUTA ELISA Kit
Reactivity	Human
Range	20-3500ng/L
Tested Applications	ELISA
Target	CUTA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	19,116
Uniprot ID	O60888
Solubility (25°C)	4000ng/L
Sensitivity	9.18ng/L
Entrez	51596
Expiration Date	Please enquire.