

Product Datasheet

Human Protein Cramped-like ELISA Kit (orb1209533)

Description	Human Protein Cramped-like, CRAMP1L ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	CRAMP1L
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	134,718
Uniprot ID	Q96RY5
Solubility (25°C)	2400ng/L
Sensitivity	5.24ng/L
Entrez	57585
Expiration Date	Please enquire.